

## Pexmetinib (ARRY-614)

Catalog No: #S6934


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pexmetinib (ARRY-614)
Brief Description	Inhibitors
Purification	/
Target Name	p38 MAPK inhibitor
Calculated MW	556.63
Formulation	C31H33FN6O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: CC1=CC=C(C=C1)N1N=C(C=C1NC(=O)NCC1=C(OC2=CC3=C(C=C2)N(CCO)N=C3)C=CC(F)=C1)C(C)(C)C  
 |c:3,5,9,11,19,24,26,33,36,39,t:1,22|Pathways: MAPK Receptor: p38 MAPK Boiling pt: Melting pt: Solubility: Ethanol: 100 mg/mL (179.65 mM) DMSO: 100 mg/mL (179.65 mM) Water:

Note: This product is for in vitro research use only