

Picolinamide

Catalog No: #S6942

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Picolinamide
Brief Description	Inhibitors
Purification	/
Target Name	PARP inhibitor
Calculated MW	122.12
Formulation	C6H6N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C1=CC=NC(=C1)C(=O)N Pathways: DNA Damage Receptor: PARP Boling pt: Melting pt: Solubility: Ethanol: 24 mg/mL (196.52 mM)

DMSO: 24 mg/mL (196.52 mM)

Water: 24 mg/mL (196.52 mM)

Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only