

Picropodophyllin (PPP)

Catalog No: #S6943


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Picropodophyllin (PPP)
Brief Description	Inhibitors
Purification	98.00%
Target Name	IGF-1R inhibitor
Calculated MW	414.41
Formulation	C22H22O8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: [H][C@]12COC(=O)[C@@]1([H])[C@H](C1=CC(OC)=C(OC)C(OC)=C1)C1=CC3=C(OCOC3)C=C1[C@@H]2O
 |r,c:20,31,t:10,14,23,25|Pathways: Tyrosine Kinase/Adaptors Receptor: IGF-1R Boling pt: Melting pt: Solubility: Ethanol: 1 mg/mL warmed (2.41 mM)
 DMSO: 82 mg/mL (197.87 mM) warming
 Water:

References

1. Lemaire M, Fristedt C, Agarwal P, et al. Clin Cancer Res. 2012 Apr 15;18(8):2230-9.

Note: This product is for in vitro research use only