


Description

Product Name	PLX7904
Brief Description	Inhibitors
Purification	/
Target Name	Raf inhibitor
Calculated MW	512.53
Formulation	C24H22F2N6O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CCN(C)S(=O)(=O)NC1=CC=C(F)C(C(=O)C2=CNC3=NC=C(C=C23)C2=CN=C(N=C2)C2CC2)=C1F
 |c:21,29,31,37,t:8,10,16,19,23,27|Pathways: MAPK Receptor: B-Raf (V600E) Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only