

NF-YC Polyclonal Antibody

Catalog No: #41223

Package Size: #41223-1 50ul #41223-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NF-YC Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	NF-YC Polyclonal Antibody detects endogenous levels of NF-YC protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human NF-YC.
Target Name	NF-YC
Other Names	NFYC; Nuclear transcription factor Y subunit gamma; CAAT box DNA-binding protein subunit C; Nuclear transcription factor Y subunit C; NF-YC; Transactivator HSM-1/2
Accession No.	Swiss-Prot: Q13952NCBI Gene ID: 4802
Uniprot	Q13952
GeneID	4802;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using NF-YC Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using NF-YC Polyclonal Antibody

Note: This product is for in vitro research use only