

Prasugrel Hydrochloride

Catalog No: #S6952


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prasugrel Hydrochloride
Brief Description	Inhibitors
Purification	/
Target Name	P2 Receptor inhibitor
Calculated MW	409.9
Formulation	C ₂₀ H ₂₀ FNO ₃ SHCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:Cl.CC(=O)OC1=CC2=C(CCN(C2)C(C(=O)C2CC2)C2=C(F)C=CC=C2)S1

|c:20,23,25,t:4,6|Pathways:NeuroscienceReceptor:P2Y12Boling pt:Melting pt:Solubility:Ethanol: 9 mg/mL (21.95 mM)

DMSO: 46 mg/mL (112.22 mM)

Water:

Note: This product is for in vitro research use only