

## Prasugrel Hydrochloride

Catalog No: #S6952


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Prasugrel Hydrochloride
Brief Description	Inhibitors
Purification	/
Target Name	P2 Receptor inhibitor
Calculated MW	409.9
Formulation	C <sub>20</sub> H <sub>20</sub> FNO <sub>3</sub> SHCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: Cl.CC(=O)OC1=CC2=C(CCN(C2)C(C(=O)C2CC2)C2=C(F)C=CC=C2)S1

|c:20,23,25,t:4,6|Pathways: Neuroscience Receptor: P2Y<sub>12</sub> Boling pt: Melting pt: Solubility: Ethanol: 9 mg/mL (21.95 mM)

DMSO: 46 mg/mL (112.22 mM)

Water:

Note: This product is for in vitro research use only