

Prilocaine

Catalog No: #S6953


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prilocaine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Sodium Channel inhibitor
Calculated MW	220.31
Formulation	C13H20N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CCCNC(C)C(=O)NC1=C(C)C=CC=C1 | c:9,12,14 | Pathways: Membrane transporter/Ion channel Receptor: Sodium Channel
Boiling pt: Melting pt: Solubility: Ethanol: 44 mg/mL (199.71 mM)
DMSO: 44 mg/mL (199.71 mM)

Water:

References

1. Eriksson E. Prilocaine. Acta Chir Scand Suppl. 1966;358:1-82.

Note: This product is for in vitro research use only