

## Propafenone HCl


Catalog No: #S6955

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Propafenone HCl
Brief Description	Inhibitors
Purification	/
Target Name	Sodium Channel inhibitor
Calculated MW	377.9
Formulation	C <sub>21</sub> H <sub>27</sub> NO <sub>3</sub> HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: CCCNCC(O)COC1=CC=CC=C1C(=O)CCC1=CC=CC=C1 | c:11,13,22,24,t:9,20 | Pathways: Membrane transporter/Ion channel Receptor: Sodium Channel Boiling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only