

## Protirelin

Catalog No: #S6956


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Protirelin
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	362.38
Formulation	C16H22N6O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Others SMILES:NC(=O)[C@@H]1CCCN1C(=O)[C@H](CC1=CNC=N1)NC(=O)[C@@H]1CCC(=O)N1 |r,c:16,t:13|Pathways:Others

Receptor:BoIing pt:Melting pt:Solubility:Ethanol: 72 mg/mL (198.68 mM)

DMSO: 72 mg/mL (198.68 mM)

Water: 72 mg/mL (198.68 mM)

Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only