


## Description

|                   |   |
|-------------------|---|
| Product Name      | PRT-060318                                      |
| Brief Description | Inhibitors                                      |
| Purification      | /   |
| Target Name       | Syk inhibitor                                   |
| Calculated MW     | 449.81  |
| Formulation       | C18H24N6O3HCl                                   |
| Storage           | 3 years -20°C powder; 2 years -80°C in solvent; |

## Images



## Product Description

Research Area: Others SMILES: CC1=CC(=CC=C1)NC2=NC(=NC=C2C(=O)N)N[C@@H]3CCCC[C@@H]3N Pathways: Angiogenesis

Receptor: Syk Boiling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only