

NFκB-p100 Polyclonal Antibody

Catalog No: #41224

Package Size: #41224-1 50ul #41224-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | NFκB-p100 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications | WB IHC ELISA |
| Species Reactivity | Hu Ms Rt |
| Specificity | NFκB-p100 Polyclonal Antibody detects endogenous levels of NFκB-p100 protein. |
| Immunogen Type | peptide |
| Immunogen Description | Synthesized peptide derived from human NFκB-p100 around the non-phosphorylation site of S869. |
| Target Name | NFκB-p100 |
| Other Names | NFKB2; LYT10; Nuclear factor NF-kappa-B p100 subunit; DNA-binding factor KBF2; H2TF1; Lymphocyte translocation chromosome 10 protein; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2; Oncogene Lyt-10; Lyt10 |
| Accession No. | Swiss-Prot: Q00653NCBI Gene ID: 4791 |
| Uniprot | Q00653 |
| GeneID | 4791; |
| SDS-PAGE MW | 97kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using NFκB-p100 Polyclonal Antibody



Western Blot analysis of Jurkat cells using NFκB-p100 Polyclonal Antibody

Note: This product is for in vitro research use only