

## Reboxetine mesylate

Catalog No: #S6963


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Reboxetine mesylate
Brief Description	Inhibitors
Purification	99.72%
Target Name	Norepinephrine inhibitor
Calculated MW	409.5
Formulation	C <sub>19</sub> H <sub>23</sub> NO <sub>3</sub> CH <sub>4</sub> O <sub>3</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:CS(O)(=O)=O.[H][C@@]1(CNCCO1)[C@H](OC1=C(OCC)C=CC=C1)C1=CC=CC=C1  
|r,c:14,19,21,26,28,t:24|Pathways:NeuroscienceReceptor:norepinephrine reuptakeBoling pt:Melting pt:Solubility:Ethanol:

## References

1. Wong EH, et al. Biol Psychiatry, 2000, 47(9), 818-829.

Note: This product is for in vitro research use only