


Description

Product Name	Scriptaid
Brief Description	Inhibitors
Purification	98.00%
Target Name	HDAC inhibitor
Calculated MW	326.35
Formulation	C18H18N2O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:ONC(=O)CCCCCN1C(=O)C2=CC=CC3=CC=CC(C1=O)=C23 |c:14,18,t:12,16,23|Pathways:Chromatin/Epigenetic;
DNA Damage/DNA Repair; NF-kbReceptor:HDACBoiling pt:Melting pt:Solubility:Ethanol:

References

1. Su GH, et al. Cancer Res, 2000, 60(12), 3137-3142.

Note: This product is for in vitro research use only