


Description

Product Name	SGI-7079
Brief Description	Inhibitors
Purification	98.00%
Target Name	c-Met/HGFR inhibitor; c-RET inhibitor; FLT inhibitor; Src inhibitor; TAM Receptor inhibitor
Calculated MW	455.53
Formulation	C26H26FN7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CN1CCN(CC1)C1=CC=C(NC2=NC3=C(C(C)=CN3)C(=N2)C2=CC(CC#N)=CC=C2)C=C1F

|c:18,22,30,32,35,t:8,10,13,15,25|Pathways: Angiogenesis; Apoptosis; Tyrosine Kinase/Adaptors/Receptor: Met; RET; FLT3; Yes; AXL; Mer/Boling

pt: Melting pt: Solubility: Ethanol:

References

1. Myers SH, et al. J Med Chem. 2016, 59(8):3593-608.

Note: This product is for in vitro research use only