

Sildenafil Mesylate

Catalog No: #S6983


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sildenafil Mesylate
Brief Description	Inhibitors
Purification	/
Target Name	PDE inhibitor
Formulation	C22H30N6O4SxCH4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CS(O)(=O)=O.CCCC1=NN(C)C2=C1N=C(NC2=O)C1=CC(=CC=C1OCC)S(=O)(=O)N1CCN(C)CC1

|c:11,14,22,24,t:7,20|Pathways: Metabolism

Receptor: Phosphodiesterase 5 Boiling pt: Melting pt: Solubility: Ethanol: 4 mg/mL DMSO: 100 mg/mL 100 Water: 100 mg/mL

100 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only