

SirReal2


Catalog No: #S6984

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SirReal2
Brief Description	Inhibitors
Purification	/
Target Name	Sirtuin inhibitor
Calculated MW	420.55
Formulation	C22H20N4OS2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC1=CC(C)=NC(SCC(=O)NC2=NC=C(CC3=CC=CC4=CC=CC=C34)S2)=N1

|c:4,19,23,30,t:1,12,14,17,21,25|Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair Receptor: SIRT2 Boling pt: Melting pt: Solubility: Ethanol:

Note: This product is for in vitro research use only