

## Sivelestat

Catalog No: #S6986


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sivelestat
Brief Description	Inhibitors
Purification	98.00%
Target Name	Serine Protease inhibitor
Calculated MW	434.46
Formulation	C20H22N2O7S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: injury SMILES: CC(C)(C)C(=O)OC1=CC=C(C=C1)S(=O)(=O)NC1=CC=CC=C1C(=O)NCC(O)=O

[c:9,11,20,22,t:7,18] Pathways: Proteases/Proteasome Receptor: neutrophil elastase Boiling pt: Melting pt: Solubility: Ethanol: 86 mg/mL (197.94 mM) DMSO: 86 mg/mL (197.94 mM) Water:

## References

1. Li G, et al. Int J Mol Med. 2016, 38(3):767-75.

Note: This product is for in vitro research use only