

Sivelestat

Catalog No: #S6986


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sivelestat
Brief Description	Inhibitors
Purification	98.00%
Target Name	Serine Protease inhibitor
Calculated MW	434.46
Formulation	C20H22N2O7S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:injurySMILES:CC(C)(C)C(=O)OC1=CC=C(C=C1)S(=O)(=O)NC1=CC=CC=C1C(=O)NCC(O)=O

[c:9,11,20,22,t:7,18]Pathways:Proteases/ProteasomeReceptor:neutrophil elastaseBoiling pt:Melting pt:Solubility:Ethanol: 86 mg/mL (197.94 mM)DMSO: 86 mg/mL (197.94 mM)Water:

References

1.Li G, et al. Int J Mol Med. 2016, 38(3):767-75.

Note: This product is for in vitro research use only