

Sodium Nitroprusside Dihydrate

Catalog No: #S6991


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium Nitroprusside Dihydrate
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	297.95
Formulation	C ₅ H ₄ FeN ₆ Na ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: [C-]#N.[C-]#N.[C-]#N.[C-]#N.[C-]#N.[N-]=O.O.O.[Na+].[Na+].[Fe+4] Pathways: Others Receptor: guanylate cyclase Boiling pt: Melting pt: Solubility: Ethanol

Note: This product is for in vitro research use only