

Sodium Tauroursodeoxycholate (TUDC)

Catalog No: #S6993


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium Tauroursodeoxycholate (TUDC)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	521.69
Formulation	C ₂₆ H ₄₄ NNaO ₆ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: [H][C@@]1(CC[C@@]2([H])[C@]3([H])[C@@H](O)C[C@]4([H])C[C@H](O)CC[C@]4(C)[C@@]3([H])CC[C@]12C)[C@H](C)CC(=O)NCCS(=O)(=O)O[Na] |r| Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol 71 mg/mL warmed (136.09 mM)

DMSO: 100 mg/mL (191.68 mM)

Water: 100 mg/mL (191.68 mM)

Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only