


Description

Product Name	STA-21
Brief Description	Inhibitors
Purification	/
Target Name	STAT inhibitor
Calculated MW	306.31
Formulation	C19H14O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC1CC2=C(C(=O)C1)C3=C(C=C2)C(=O)C4=C(C3=O)C=CC=C4OPathways:JAK/STAT
;Stem Cells & Wnt

Receptor:STAT3Boiling pt:Melting pt:Solubility:DMSO: Water: Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only