

Tedizolid Phosphate

Catalog No: #S6998


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tedizolid Phosphate
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	450.32
Formulation	C17H16FN6O6P
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CN1N=NC(=N1)C1=NC=C(C=C1)C1=CC=C(C=C1F)N1C[C@H](COP(O)(O)=O)OC1=O
|r,c:2,4,9,11,16,18,t:7,14|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol

Note: This product is for in vitro research use only