

Tiagabine hydrochloride


Catalog No: #S7000

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tiagabine hydrochloride
Brief Description	Inhibitors
Purification	/
Target Name	GABA Receptor inhibitor
Calculated MW	412.01
Formulation	C ₂₀ H ₂₅ NO ₂ S ₂ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: Cl.CC1=C(SC=C1)C(=CCCN1CCC[C@H](C1)C(O)=O)C1=C(C)C=CS1 |r,c:4,21,24,t:1|Pathways: Membrane transporter/Ion channel; Neuroscience Receptor: GABA uptake Boling pt: Melting pt: Solubility: Ethanol 82 mg/mL (199.02 mM) DMSO: 82 mg/mL (199.02 mM) Water: 11 mg/mL (26.69 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only