

TIC10


Catalog No: #S7001

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TIC10
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	386.49
Formulation	C24H26N4O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC1=CC=CC=C1CN1C2=NCCN2C2=C(CN(CC3=CC=CC=C3)CC2)C1=O
|c:3,5,23,25,t:1,10,16,21|Pathways:OthersReceptor:TRAILBoling pt:Melting pt:Solubility:Ethanol

Note: This product is for in vitro research use only