

Trihexyphenidyl hydrochloride

Catalog No: #S7002


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trihexyphenidyl hydrochloride
Brief Description	Inhibitors
Purification	/
Target Name	AChR inhibitor
Calculated MW	337.93
Formulation	C ₂₀ H ₃₁ NOHCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: Cl.Oc(CCN1CCCCC1)(C1CCCCC1)C1=CC=CC=C1 |c:20,22,t:18|Pathways: Neuroscience Receptor: AChR Boling
 pt: Melting pt: Solubility: Ethanol 24 mg/mL warming (71.02 mM) DMSO: 7 mg/mL warming (20.71 mM) Water:

Note: This product is for in vitro research use only