

## NFκB-p65 Polyclonal Antibody

Catalog No: #41228

Package Size: #41228-1 50ul #41228-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NFκB-p65 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	NFκB-p65 Polyclonal Antibody detects endogenous levels of NFκB-p65 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human NFκB-p65 around the non-phosphorylation site of S536.
Target Name	NFκB-p65
Other Names	RELA; NFKB3; Transcription factor p65; Nuclear factor NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3
Accession No.	Swiss-Prot: Q04206NCBI Gene ID: 5970
Uniprot	Q04206
GeneID	5970;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of Jurkat cells using NFκB-p65 Polyclonal Antibody



Western Blot analysis of COLO 293T HELA cells using NFκB-p65 Polyclonal Antibody

---

Note: This product is for in vitro research use only