

## Vacquinol-1

Catalog No: #S7008


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Vacquinol-1
Brief Description	Inhibitors
Purification	98.00%
Target Name	JNK agonist
Calculated MW	425.78
Formulation	C <sub>21</sub> H <sub>21</sub> ClN <sub>2</sub> O · 2HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: C1CCNC(C1)C(C2=CC(=NC3=CC=CC=C32)C4=CC=C(C=C4)Cl)O Pathways: MAPK Receptor: MKK4 Boling pt: Melting pt: Solubility: Ethanol

## References

1. Sander P., Vacquinol-1 inducible cell death in glioblastoma multiforme is counter regulated by TRPM7 activity induced by exogenous ATP. *Oncotarget*. 2017 May 23;8(21):35124-35137.

Note: This product is for in vitro research use only