


## Description

Product Name	VO-Ohpic
Brief Description	Inhibitors
Purification	99.00%
Target Name	PTEN inhibitor
Calculated MW	415.2
Formulation	C12H9N2O8V3H2OH
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: C1=CC(=C(N=C1)C(=O)O)O.C1=CC(=C(N=C1)C(=O)O)O.O.O.O.O=[VH] Pathways: Others Receptor: PTEN Boling  
 pt: Melting pt: Solubility: Ethanol

Note: This product is for in vitro research use only