

VPS34-IN1

Catalog No: #S7015


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	VPS34-IN1
Brief Description	Inhibitors
Purification	/
Target Name	PI3K inhibitor
Calculated MW	425.91
Formulation	C21H24ClN7O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC(C)(O)CNC1=NC=C(C(CC2CC2)=N1)C1=CC=NC(NC2=CC=NC(Cl)=C2)=N1

|c:8,15,20,26,29,31,t:6,18,24|Pathways: PI3K/Akt/mTOR signaling Receptor: Vps34 Boling pt: Melting pt: Solubility: Ethanol 85 mg/mL (199.57 mM) DMSO: 85 mg/mL (199.57 mM) Water:

Note: This product is for in vitro research use only