

ANTHRAQUINONE-2-CARBOXYLIC ACID

Catalog No: #S7035

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ANTHRAQUINONE-2-CARBOXYLIC ACID
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antiviral; Syk
Calculated MW	252.22
Formulation	C15H8O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Others SMILES: C1=CC=C2C(=C1)C(=O)C3=C(C2=O)C=C(C=C3)C(=O)O

Pathways: Angiogenesis; Immunology/Inflammation Receptor: IRAK-1; Syk Boiling pt: Melting pt: Solubility: N/A Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Park JG, et al.

Anti-Inflammatory and Antinociceptive Activities of Anthraquinone-2-Carboxylic Acid. *Mediators Inflamm.* 2016;2016:1903849.

Note: This product is for in vitro research use only