

4-METHYL-6,7-METHYLENEDIOXYCOUMARIN

Catalog No: #S7038

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-METHYL-6,7-METHYLENEDIOXYCOUMARIN
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	204.18
Formulation	C11H8O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:OthersSMILES:CC1=CC(=O)OC2=CC3=C(C=C12)OCO3Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:N/AAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only