

Angiotensin II human

Catalog No: #S7040

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin II human
Brief Description	Inhibitors
Purification	98.00%
Target Name	RAAS
Calculated MW	1046.18
Formulation	C50H71N13O12
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research

Area: cancer
 SMILES: CC[C@H](C)[C@@H](C(=O)N[C@@H](CC1=CN=CN1)C(=O)N2CCC[C@H]2C(=O)N[C@@H](CC3=CC=CC=C3)C(=O)O)NC(=O)[C@H](CC4=CC=C(C=C4)O)NC(=O)[C@H](C(C)C)NC(=O)[C@H](CCCN=C(N)N)NC(=O)[C@H](CC(=O)O)N
 Pathways: Endocrinology/Hormones
 Receptor: AT receptor
 Boiling pt: Melting pt: Solubility: H₂O: 50 mg/mL
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. de Gasparo M, et al. International union of pharmacology. XXIII. The angiotensin II receptors. *Pharmacol Rev.* 2000 Sep;52(3):415-72.

Note: This product is for in vitro research use only