

NTE Polyclonal Antibody

Catalog No: #41256

Package Size: #41256-1 50ul #41256-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NTE Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	NTE Polyclonal Antibody detects endogenous levels of NTE protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human NTE.
Target Name	NTE
Other Names	PNPLA6; NTE; Neuropathy target esterase; Patatin-like phospholipase domain-containing protein 6
Accession No.	Swiss-Prot: Q8IY17NCBI Gene ID: 10908
Uniprot	Q8IY17
GeneID	10908;
SDS-PAGE MW	150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of SH-SY5Y cells using NTE Polyclonal Antibody

Note: This product is for in vitro research use only