

NUB1 Polyclonal Antibody

Catalog No: #41257

Package Size: #41257-1 50ul #41257-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NUB1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Mk
Specificity	NUB1 Polyclonal Antibody detects endogenous levels of NUB1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human NUB1.
Target Name	NUB1
Other Names	NUB1; NYREN18; NEDD8 ultimate buster 1; Negative regulator of ubiquitin-like proteins 1; Renal carcinoma antigen NY-REN-18
Accession No.	Swiss-Prot: Q9Y5A7NCBI Gene ID: 51667
Uniprot	Q9Y5A7
GeneID	51667;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of COS7 cells using NUB1 Polyclonal Antibody



Note: This product is for in vitro research use only