

## Olfactory receptor 10T2 Polyclonal Antibody

Catalog No: #41265

Package Size: #41265-1 50ul #41265-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Olfactory receptor 10T2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 10T2 Polyclonal Antibody detects endogenous levels of Olfactory receptor 10T2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 10T2.
Target Name	Olfactory receptor 10T2
Other Names	OR10T2; Olfactory receptor 10T2; Olfactory receptor OR1-3
Accession No.	Swiss-Prot: Q8NGX3NCBI Gene ID: 128360
Uniprot	Q8NGX3
GeneID	128360;
SDS-PAGE MW	30kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293T-UV cells using Olfactory receptor 10T2 Polyclonal Antibody

---

Note: This product is for in vitro research use only