

Olfactory receptor 11G2 Polyclonal Antibody

Catalog No: #41266

Package Size: #41266-1 50ul #41266-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 11G2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 11G2 Polyclonal Antibody detects endogenous levels of Olfactory receptor 11G2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 11G2.
Target Name	Olfactory receptor 11G2
Other Names	OR11G2; Olfactory receptor 11G2; Olfactory receptor OR14-34
Accession No.	Swiss-Prot: Q8NGC1NCBI Gene ID: 390439
Uniprot	Q8NGC1
GeneID	390439;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using Olfactory receptor 11G2 Polyclonal Antibody

Note: This product is for in vitro research use only