

## Olfactory receptor 4C16 Polyclonal Antibody

Catalog No: #41271

Package Size: #41271-1 50ul #41271-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Olfactory receptor 4C16 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 4C16 Polyclonal Antibody detects endogenous levels of Olfactory receptor 4C16 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4C16.
Target Name	Olfactory receptor 4C16
Other Names	OR4C16; Olfactory receptor 4C16; Olfactory receptor OR11-135
Accession No.	Swiss-Prot: Q8NGL9NCBI Gene ID: 219428
Uniprot	Q8NGL9
GeneID	219428;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of SH-SY5Y cells using Olfactory receptor 4C16 Polyclonal Antibody

---

Note: This product is for in vitro research use only