

Olfactory receptor 4K14 Polyclonal Antibody

Catalog No: #41274

Package Size: #41274-1 50ul #41274-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 4K14 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Mk
Specificity	Olfactory receptor 4K14 Polyclonal Antibody detects endogenous levels of Olfactory receptor 4K14 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4K14.
Target Name	Olfactory receptor 4K14
Other Names	OR4K14; Olfactory receptor 4K14; Olfactory receptor OR14-22
Accession No.	Swiss-Prot: Q8NGD5NCBI Gene ID: 122740
Uniprot	Q8NGD5
GeneID	122740;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using Olfactory receptor 4K14 Polyclonal Antibody

Note: This product is for in vitro research use only