

Olfactory receptor 4Q3 Polyclonal Antibody

Catalog No: #41275

Package Size: #41275-1 50ul #41275-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 4Q3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 4Q3 Polyclonal Antibody detects endogenous levels of Olfactory receptor 4Q3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4Q3.
Target Name	Olfactory receptor 4Q3
Other Names	OR4Q3; C14orf13; OR4Q4; Olfactory receptor 4Q3; Olfactory receptor 4Q4; Olfactory receptor OR14-3
Accession No.	Swiss-Prot: Q8NH05NCBI Gene ID: 441669
Uniprot	Q8NH05
GeneID	441669;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using Olfactory receptor 4Q3 Polyclonal Antibody

Note: This product is for in vitro research use only