

## Olfactory receptor 51B2 Polyclonal Antibody

Catalog No: #41276

Package Size: #41276-1 50ul #41276-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Olfactory receptor 51B2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 51B2 Polyclonal Antibody detects endogenous levels of Olfactory receptor 51B2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Olfactory receptor 51B2.
Target Name	Olfactory receptor 51B2
Other Names	OR51B2; OR51B1P; Olfactory receptor 51B2; Odorant receptor HOR5 $\beta$ 3; Olfactory receptor 51B1
Accession No.	Swiss-Prot: Q9Y5P1NCBI Gene ID: 79345
Uniprot	Q9Y5P1
GeneID	79345;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of HT29 cells using Olfactory receptor 51B2 Polyclonal Antibody

---

Note: This product is for in vitro research use only