

## Olfactory receptor 52E4 Polyclonal Antibody

Catalog No: #41277

Package Size: #41277-1 50ul #41277-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Olfactory receptor 52E4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 52E4 Polyclonal Antibody detects endogenous levels of Olfactory receptor 52E4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Olfactory receptor 52E4.
Target Name	Olfactory receptor 52E4
Other Names	OR52E4; Olfactory receptor 52E4; Olfactory receptor OR11-55
Accession No.	Swiss-Prot: Q8NGH9NCBI Gene ID: 390081
Uniprot	Q8NGH9
GeneID	390081;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of HeLa cells using Olfactory receptor 52E4 Polyclonal Antibody

---

Note: This product is for in vitro research use only