

Olfactory receptor 52K1 Polyclonal Antibody

Catalog No: #41279

Package Size: #41279-1 50ul #41279-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 52K1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 52K1 Polyclonal Antibody detects endogenous levels of Olfactory receptor 52K1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Olfactory receptor 52K1.
Target Name	Olfactory receptor 52K1
Other Names	OR52K1; Olfactory receptor 52K1; Olfactory receptor OR11-8
Accession No.	Swiss-Prot: Q8NGK4NCBI Gene ID: 390036
Uniprot	Q8NGK4
GeneID	390036;
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using Olfactory receptor 52K1 Polyclonal Antibody

Note: This product is for in vitro research use only