

## Olfactory receptor 5I1 Polyclonal Antibody

Catalog No: #41281

Package Size: #41281-1 50ul #41281-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Olfactory receptor 5I1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 5I1 Polyclonal Antibody detects endogenous levels of Olfactory receptor 5I1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Olfactory receptor 5I1.
Target Name	Olfactory receptor 5I1
Other Names	OR5I1; OLF1; Olfactory receptor 5I1; Olfactory receptor OR11-159; Olfactory receptor-like protein OLF1
Accession No.	Swiss-Prot: Q13606NCBI Gene ID: 10798
Uniprot	Q13606
GeneID	10798;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293T cells using Olfactory receptor 5I1 Polyclonal Antibody

---

Note: This product is for in vitro research use only