

Olfactory receptor 5M3 Polyclonal Antibody

Catalog No: #41282

Package Size: #41282-1 50ul #41282-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 5M3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 5M3 Polyclonal Antibody detects endogenous levels of Olfactory receptor 5M3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Olfactory receptor 5M3.
Target Name	Olfactory receptor 5M3
Other Names	OR5M3; Olfactory receptor 5M3; Olfactory receptor OR11-191
Accession No.	Swiss-Prot: Q8NGP4NCBI Gene ID: 219482
Uniprot	Q8NGP4
GeneID	219482;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HEPG2 cells using Olfactory receptor 5M3 Polyclonal Antibody

Note: This product is for in vitro research use only