

## Olfactory receptor 5W2 Polyclonal Antibody

Catalog No: #41283

Package Size: #41283-1 50ul #41283-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Olfactory receptor 5W2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 5W2 Polyclonal Antibody detects endogenous levels of Olfactory receptor 5W2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Olfactory receptor 5W2.
Target Name	Olfactory receptor 5W2
Other Names	OR5W2; OR5W2P; OR5W3P; Olfactory receptor 5W2; Olfactory receptor 5W3; Olfactory receptor OR11-155
Accession No.	Swiss-Prot: Q8NH69NCBI Gene ID: 390148
Uniprot	Q8NH69
GeneID	390148;
SDS-PAGE MW	34kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293T-UV cells using Olfactory receptor 5W2 Polyclonal Antibody

---

Note: This product is for in vitro research use only