

Olfactory receptor 6Q1 Polyclonal Antibody

Catalog No: #41284

Package Size: #41284-1 50ul #41284-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Olfactory receptor 6Q1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 6Q1 Polyclonal Antibody detects endogenous levels of Olfactory receptor 6Q1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Olfactory receptor 6Q1.
Target Name	Olfactory receptor 6Q1
Other Names	OR6Q1; Olfactory receptor 6Q1; Olfactory receptor OR11-226
Accession No.	Swiss-Prot: Q8NGQ2NCBI Gene ID: 219952
Uniprot	Q8NGQ2
GeneID	219952;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using Olfactory receptor 6Q1 Polyclonal Antibody

Note: This product is for in vitro research use only