

p115 Polyclonal Antibody

Catalog No: #41294

Package Size: #41294-1 50ul #41294-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	p115 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	p115 Polyclonal Antibody detects endogenous levels of p115 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human p115.
Target Name	p115
Other Names	USO1; VDP; General vesicular transport factor p115; Protein USO1 homolog; Transcytosis-associated protein; TAP; Vesicle-docking protein
Accession No.	Swiss-Prot: O60763NCBI Gene ID: 8615
Uniprot	O60763
GeneID	8615;
SDS-PAGE MW	108kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using p115 Polyclonal Antibody



Western Blot analysis of K562 293 cells using p115 Polyclonal Antibody

Note: This product is for in vitro research use only