

## p140Cap Polyclonal Antibody

Catalog No: #41295

Package Size: #41295-1 50ul #41295-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	p140Cap Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	p140Cap Polyclonal Antibody detects endogenous levels of p140Cap protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human p140Cap.
Target Name	p140Cap
Other Names	SRCIN1; KIAA1684; P140; SRC kinase signaling inhibitor 1; SNAP-25-interacting protein; SNIP; p130Cas-associated protein; p140Cap
Accession No.	Swiss-Prot: Q9C0H9NCBI Gene ID: 80725
Uniprot	Q9C0H9
GeneID	80725;
SDS-PAGE MW	112kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of HeLa cells using p140Cap Polyclonal Antibody



Western Blot analysis of 293T HELA cells using p140Cap Polyclonal Antibody

---

Note: This product is for in vitro research use only