

P2RY8 Polyclonal Antibody

Catalog No: #41301

Package Size: #41301-1 50ul #41301-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2RY8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	P2RY8 Polyclonal Antibody detects endogenous levels of P2RY8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human P2RY8.
Target Name	P2RY8
Other Names	P2RY8; P2Y purinoceptor 8; P2Y8
Accession No.	Swiss-Prot: Q86VZ1NCBI Gene ID: 286530
Uniprot	Q86VZ1
GeneID	286530;
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using P2RY8 Polyclonal Antibody

Note: This product is for in vitro research use only